

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human TNFRSF25 / DR3 / TNFRSF12 Protein (Fc tag)(Fc tag)

产品货号: TP07935

第三版

别名:APO-3;DDR3;DR3;LARD;TNFRSF12;TR3;TRAMP;WSL-1;WSL-LR

描述:Tumor necrosis factor receptor superfamily, member 25 (TNFRSF25), also known as Death receptor 3 (DR3) or TNFRSF12 is a member of the TNF-receptor superfamily. This receptor is expressed preferentially in the tissues enriched in lymphocytes, and it may play a role in regulating lymphocyte homeostasis. TNFRSF25/DR3/TNFRSF12 has been shown to stimulate NF-kappa B activity and regulate cell apoptosis. The signal transduction of this receptor is mediated by various death domain containing adaptor proteins. Knockout studies in mice suggested the role of this gene in the removal of self-reactive T cells in the thymus. Multiple alternatively spliced transcript variants of this gene encoding distinct isoforms have been reported, most of which are potentially secreted molecules. The alternative splicing of this TNFRSF25 encoding gene in B and T cells encounters a programmed change upon T-cell activation, which predominantly produces full-length, membrane bound isoforms, and is thought to be involved in controlling lymphocyte proliferation induced by T-cell activation.

配方:PBS

分子量:46 kDa

序列:Met 1-Gln 199

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.