

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse KLK1 / Kallikrein 1 Protein (His tag)

产品货号: TP08002

第三版

别名:0610007D04Rik;Kal;Klk1b6;Klk6;mGk-6

描述:Human tissue kallikrein (KLK1) is a serine protease, component of the KKS that has been demonstrated to exert pleiotropic beneficial effects in protection from tissue injury through its antiinflammatory, antiapoptotic, antifibrotic and antioxidative actions. Polymorphism of the human tissue kallikrein 1 (KLK1) A1789G gene is associated with susceptibility to hypertension.

配方:PBS

分子量:28.3 kDa

序列:Met1-Asp261

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.