

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Rat Thrombomodulin / THBD Protein (His tag)

产品货号: TP08038

第三版

别名:THBD

描述:Thrombomodulin, also known as THBD(CD141), is an integral membrane protein which reduces blood coagulation by converting thrombin to an anticoagulant enzyme from a procoagulant enzyme. Thrombomodulin is expressed on the surface of endothelial cells and serves as a cofactor for thrombin. It is also expressed on human mesothelial cell, monocyte and a dendritic cell subset. Thrombomodulin functions as a cofactor in the thrombin-induced activation of protein C in the anticoagulant pathway by forming a 1:1 stoichiometric complex with thrombin. Thrombomodulin also regulates C3b inactivation by factor I. Mutations in the thrombomodulin gene have also been reported to be associated with atypical hemolytic-uremic syndrome.

配方:PBS

分子量:55.2 kDa

序列:Met1-Ser517

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.