
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human SMOC1 Protein (His tag)

产品货号: TP08042

第三版

别名:OAS

描述:SPARC-related modular calcium-binding protein 1, also known as secreted modular calcium-binding protein 1 and SMOC1, is a member of the SPARC family. SMOC1 is widely expressed in many tissues with a strongest signal in ovary. It contains two EF-hand domains, one Kazal-like domain and two thyroglobulin type-1 domains. Extracellular matrix proteins have been implicated in the regulation of osteoblast differentiation of bone marrow derived mesenchymal stem cells (BMSCs) through paracrine or autocrine mechanisms. SMOC1 is a regulator of osteoblast differentiation of BMSCs. SMOC1 is highly expressed and secreted in BMSCs stimulated with osteogenic medium (OSM). SMOC1 and SMOC2 are matricellular proteins thought to influence growth factor signaling, migration, proliferation, and angiogenesis. SMOC1 and SMOC2 may mediate intercellular signaling and cell type-specific differentiation during gonad and reproductive tract development.

配方:PBS

分子量:47 kDa

序列:Met 1-Val 435

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.