

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Granzyme B / GZMB Protein (His tag)

产品货号: TP08131

第三版

别名:CCPI;CGL-1;CGL1;CSP-B;CSPB;CTLA1;CTSG1;HLP;SECT

**描述:**Granzyme B, also known as GZMB, is the most prominent member of the granzyme family of cell death-inducing serine proteases expressed in the granules of cytotoxic T lymphocytes (CTLs) and NK cells. Granzyme B enters the target cells depending on another membrane-binding granule protein, perforin, results in the activation of effector caspases and mitochondrial depolarization through caspase-dependent and -independent pathways, and consequently induces rapid cell apoptosis. Over 30 substrates of GZMB have been identified including the key substrate caspase-3, ICAD and Bid. GZMB is suggested to protect the host by lysing cells bearing on their surface 'nonself' antigens such as bacterial and viral infected-cells and tumor cells, and accordingly plays an essential role in immunosurveillance.

**配方:**PBS

**分子量:**27 kDa

**序列:**Met 1-Tyr 247

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.