



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human DCBLD2 / ESDN / CLCP1 Protein (His tag)

产品货号: TP08141

第三版

别名:CLCP1;ESDN

描述:DCBLD2, also known as ESDN and CLCP1, localizes in various compartments. DCBLD2 is up-regulated in lung cancers and is regulated by transcription factor AP-2 alpha (TFAP2A), a component of activator protein-2 (AP-2) that is known to regulate IL-8 production in human lung fibroblasts and epithelial cells. DCBLD2 could be related to FEV(1)-related phenotypes in asthmatics. DCBLD2 gene is expressed at very high level. DCBLD2 is proposed to participate in processes such as intracellular receptor mediated signaling pathway, negative regulation of cell growth and so on.

配方:PBS

分子量:52.6 kDa

序列:Met 1-Ala482

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.