

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human CD32b / FCGR2B Protein (His & AVI tag)

产品货号: TP08233

第三版

别名:CD32;CD32B;Fc gamma RIIB;FCGR2

描述:MD-1 and MD-2 are secretory glycoproteins that exist on the cell surface in complexes with transmembrane proteins. MD-1 is anchored by radioprotective 105 (RP105) which is a molecule containing leucine-rich repeats and is expressed on B cells, dendritic cells and macrophages, while MD-2 is associated with TLR4. MD-1 is required for efficient RP105 cell surface expression and function. It is indicated that the RP105/MD1 complex, in conjunction with TLR4, mediates the innate immune response to LPS in B cells, and also plays a role in protecting against apoptosis, B-cell proliferation, etc. Mouse MD-1 cDNA encodes a 162 amino acid precursor protein with a putative 19 aa signal peptide and two potential N-linked glycosylation sites. It shares 40% and 66% amino acid sequence identity with chicken and human MD-1 respectively. MD-1 is mainly expressed in spleen, and also detectable in liver, brain, thymus, and kidney.

配方:PBS

分子量:24 kDa

序列:Ala 46-Ile 224

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.