

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human DBH / Dopamine beta-Hydroxylase Protein (His tag)

产品货号: TP08269

第三版

别名:DBM

描述:DBH is a 290 kDa copper-containing oxygenase. It can be detected in noradrenergic nerve terminals of the central and peripheral nervous systems, and is also expressed in chromaffin cells of the adrenal medulla. DBH contains our identical subunits, and its activity requires ascorbate as a cofactor. It functions in the synthesis of small-molecule neurotransmitters that is membrane-bound, making norepinephrine the only transmitter synthesized inside vesicles. DBH has been shown to be associated with decision making and addictive behaviors such as alcohol and smoking, attention deficit hyperactivity disorder, and also with neurological diseases such as Schizophrenia and Alzheimer's.

配方:PBS

分子量:67.3 kDa

序列:Ser26-Gly603

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储:+4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.