
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse HAI-2 / SPINT2 Protein (His tag)

产品货号: TP08272

第三版

别名:AL024025;C76321;HAI-2

描述:Choanal (CA) and gastrointestinal atresias (GA) are an important feature of syndromic congenital sodium diarrhea (sCSD), a disorder recently associated with mutations in the gene for serine protease inhibitor type 2 (SPINT2). The SPINT2 gene is epigenetically silenced or downregulated in human cancers, altering the balance of HGF activation/inhibition ratio, which contributes to cancer development and progression. SPINT2 is a tumor suppressor gene that inhibits proteases implicated in cancer progression, like HGFA, hepsin and matriptase. Loss of SPINT2 expression in tumors has been associated with gene promoter hypermethylation. SPINT2 (serine peptidase inhibitor Kunitz type 2), a proteolytic inhibitor of hepatocyte growth factor activator (HGFA), which has a significant role in the suppression of the HGF-MET pathway and malignant melanoma progression.

配方:PBS

分子量:14 kDa

序列:Met 1-Lys 140

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.