
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human SPHK1 / Sphingosine Kinase 1 Protein (His & GST tag)

产品货号: TP08332

第三版

别名:SPHK;SPHK1;SPK

描述:SPHK1, also known as sphingosine Kinase 1, catalyzes the phosphorylation of sphingosine to form sphingosine-1-phosphate (S1P). S1P exhibits a broad spectrum of biological activities including cell proliferation, survival, migration, cytoskeletal organization, and morphogenesis. It is a ligand for cell surface G protein-coupled receptors. SPHK 1 is a potential therapeutic target for the control of cancer and inflammation. SPHK1 plays a key role in TNF-alpha signaling and the NF-kappa-B activation pathway important in inflammatory, antiapoptotic, and immune processes.

配方:PBS

分子量:70.3 kDa

序列:Met 1-Leu 384

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.