



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human CMPK1 Protein (His tag)

产品货号: TP08338

第三版

别名:CK;CMK;CMPK;UMK;UMP-CMPK;UMPK

描述:CMPK1 plays a key role in the maintenance of pyrimidine nucleotide pool profile and for the metabolism of pyrimidine analogs in cells. It catalyzes the phosphoryl transfer from ATP to UMP, CMP, and deoxy-CMP (dCMP), resulting in the formation of ADP and the corresponding nucleoside diphosphate. CMPK1 also has a significant role in the activation of pyrimidine analogs, which are clinically useful anti-cancer and anti-viral drugs. In the meanwhile, CMPK1 functions in cellular nucleic acid biosynthesis.

配方:PBS

分子量:24.5 kDa

序列:Met 1-Gly 196

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.