

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human STK24 / MST3 Protein (His tag)

产品货号: TP08348

第三版

别名:HEL-S-95;MST3;MST3B;RP11-111L24.5;STE20;STK3

描述:MST3 encodes a serine/threonine protein kinase that functions upstream of mitogen-activated protein kinase (MAPK) signaling. The encoded protein is cleaved into two chains by caspases; the N-terminal fragment (MST3/N) translocates to the nucleus and promotes programmed cells death. There is a pseudogene for this gene on chromosome X. Alternative splicing results in multiple transcript variants.

配方:PBS

分子量:50.3 kDa

序列:Met 1-His 431

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.