
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human ADK Protein (His & GST tag)

产品货号: TP08349

第三版

别名:AK

描述:Adenosine kinase(ADK) belongs to the family of transferases. Adenosine kinase (ADK) is the key enzyme in adenosine metabolism and catalyzes ATP and adenosine into two products: ADP and AMP. Two isoforms of the enzyme adenosine kinase (ADK), which differ at their N-terminal ends, are found in mammalian cells. It has been shown that the two ADK isoforms differ only in their first exons and the promoter regions; hence they arise via differential splicing of their first exons with the other exons common to both isoforms. In adult brain, ADK is primarily present in astrocytes. Several lines of experimental evidence support a critical role of ADK in different types of brain injury associated with astrogliosis, which is also a prominent morphologic feature of temporal lobe epilepsy (TLE). It has been suggested that dysregulation of ADK in astrocytes is a common pathologic hallmark of TLE. Moreover, in vitro data suggest the existence of an additional layer of modulatory crosstalk between the astrocyte-based adenosine cycle and inflammation. ADK also contributes to CK homeostasis in vivo.

配方:PBS

分子量:68 kDa

序列:Met 1-His 345

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.