
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Casein Kinase 1 gamma 2 Protein (His Tag)

产品货号: TP08355

第三版

别名:CK1g2

描述:Casein kinase I gamma 2 isoform (CSNK1G2), a member of the large casein kinase I (CKI) subfamily, protein kinase superfamily. It may affect the development of brain, and associate with vesicular trafficking and neurotransmitter releasing from small synaptic vesicles. The CKI family includes several other isoforms (alpha, beta, gamma, and delta). Dishevelled (Dsh), another positive component of the Wnt pathway, becomes phosphorylated in response to Wnt signals. All the CKI isoforms, with the exception of gamma, increase the phosphorylation of Dsh in vivo. Casein kinase 1 gamma (CK1gamma, or CSNK1G) is associated with the cell membrane and binds to LRP. CK1gamma was found to be needed for Wnt signaling through Wnt receptor LRP. CSNK1G2 inhibits Smad3-mediated TGF-beta responses including induction of target genes and cell growth arrest, and this inhibition is dependent on CSNK1G2 kinase activity. The overexpression of CSNK1G2 in human cancers, may act as an oncoprotein during tumorigenesis. In addition, as an MTA1s-binding protein, CSNK1G2 could further potentiate the estrogen receptor (ER) corepressive function of MTA1s.

配方:PBS

分子量:47.8 kDa

序列:Met 18-Lys 415

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.