



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human TLK1 / PKU-beta Protein (His & GST tag)

产品货号: TP08357

第三版

别名:KIAA0137;PKU-beta

描述:The protein encoded by this gene is a serine/threonine kinase that may be involved in the regulation of chromatin assembly. The encoded protein is only active when it is phosphorylated, and this phosphorylation is cell cycle-dependent, with the maximal activity of this protein coming during S phase. The catalytic activity of this protein is diminished by DNA damage and by blockage of DNA replication. Three transcript variants encoding different isoforms have been found for this gene.

配方:PBS

分子量:114 kDa

序列:Met 1-Tyr 766

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.