
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human MAP4K5 / MEK5 Protein (His & GST tag)

产品货号: TP08374

第三版

别名:GCKR;KHS;KHS1;MAPKKKK5

描述:Mitogen-activated protein kinase kinase kinase 5, also known as Kinase homologous to SPS1/STE20, MAPK/ERK kinase kinase 5, MEK kinase kinase 5 and MAP4K5, is a cytoplasm protein which belongs to the protein kinase superfamily, STE Ser/Thr protein kinase family and STE20 subfamily. MAP4K5 is ubiquitously expressed in all tissues examined, with high levels in the ovary, testis and prostate. It contains one CNH domain and one protein kinase domain. MAP4K5 is highly similar to yeast SPS1/STE20 kinase. Yeast SPS1/STE20 functions near the beginning of the MAP kinase signal cascades that is essential for yeast pheromone response. MAP4K5 has been shown to interact with CRKL and TRAF2. This kinase was shown to activate Jun kinase in mammalian cells. MAP4K5 is an early component of MAP kinase signal cascades. It may play a role in the response to environmental stress. MAP4K5 appears to act upstream of the JUN N-terminal pathway.

配方:PBS

分子量:123 kDa

序列:Met 1-Tyr 846

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.