
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human TLK2 / PKU-ALPHA Protein (aa 397-772, His & GST tag)

产品货号: TP08387

第三版

别名:HsHPK;PKU-ALPHA

描述: Serine / threonine-protein kinase tousled-like 2, also known as PKU-alpha, Tousled-like kinase 2 and TLK2, is a nucleus protein which belongs to the protein kinase superfamily and Ser/Thr protein kinase family. The tousled-like kinases are an evolutionarily conserved family of proteins implicated in DNA repair, DNA replication and mitosis in metazoans and plants. Their absence from the yeasts and other eukaryotic 'microbes' suggests a specific role for them in the development of multicellular organisms. Tousled-like kinase 2 / TLK2 is widely expressed. It is present in fetal placenta, liver, kidney, pancreas, heart and skeletal muscle. It is also found in adult cell lines. Tousled-like kinase 2 / TLK2 contains one protein kinase domain. Tousled-like kinase 2 / TLK2 is rapidly and transiently inhibited by phosphorylation following the generation of DNA double-stranded breaks during S-phase. This is cell cycle checkpoint and ATM-pathway dependent and appears to regulate processes involved in chromatin assembly.

配方: PBS

分子量: 71.3 kDa

序列: Leu 397-Asn 772

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.