
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human TAOK3 / JIK / MAP3K18 Protein (aa 1-411, His & GST tag)

产品货号: TP08389

第三版

别名:DPK;JIK;MAP3K18

描述: Serine/threonine-protein kinase TAO3, also known as cutaneous T-cell lymphoma-associated antigen HD-CL-09, CTCL-associated antigen HD-CL-09, Dendritic cell-derived protein kinase, JNK / SAPK-inhibitory kinase, Jun kinase-inhibitory kinase, Kinase from chicken homolog A, Thousand and one amino acid protein 3, JIK, TAOK3 and MAP3K18, is cytoplasm and peripheral membrane protein which belongs to the protein kinase superfamily, STE Ser/Thr protein kinase family and STE20 subfamily. Protein kinases are pivotal regulators of cell signaling that modulate each other's functions and activities through site-specific phosphorylation events. TAOK3 / JIK contains one protein kinase domain. TAOK3 / JIK is ubiquitously expressed at a low level, and highly expressed in peripheral blood leukocytes (PBLs), thymus, spleen, kidney, skeletal muscle, heart and liver. TAOK3 / JIK inhibits the basal activity of Jun kinase. It is negatively regulated by epidermal growth factor (EGF). When overexpressed, TAOK3 / JIK may activate ERK1 / ERK2 and JNK / SAPK.

配方: PBS

分子量: 74.3 kDa

序列: Met 1-Asp 411

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.