

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human STK4 / MST1 Protein (His tag)

产品货号: TP08390

第三版

别名:KRS2;MST1;TIIAC;YSK3

描述:MST1 protein encoded by this gene contains four kringle domains and a serine protease domain, similar to that found in hepatic growth factor. Despite the presence of the serine protease domain, the encoded protein may not have any proteolytic activity. The receptor for this protein is RON tyrosine kinase, which upon activation stimulates ciliary motility of ciliated epithelial lung cells. This protein is secreted and cleaved to form an alpha chain and a beta chain bridged by disulfide bonds.

配方:PBS

分子量:58 kDa

序列:Glu 2-Phe 487

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.