

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human AKT2 / protein kinase B / PKB beta Protein (His & GST tag)

产品货号: TP08407

第三版

别名:HIHGHH;PKBB;PKBBETA;PRKBB;RAC-BETA

**描述:**AKT2 (AKT Serine/Threonine Kinase 2) is a Protein Coding gene. Diseases associated with AKT2 include Hypoinsulinemic Hypoglycemia With Hemihypertrophy and Diabetes Mellitus, Noninsulin-Dependent. Among its related pathways are Constitutive Signaling by AKT1 E17K in Cancer and RANK Signaling in Osteoclasts. Gene Ontology (GO) annotations related to this gene include transferase activity, transferring phosphorus-containing groups and protein tyrosine kinase activity. An important paralog of this gene is AKT1.

**配方:**PBS

**分子量:**83.6 kDa

**序列:**Met 1-Glu481

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.