
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human RSK4 / RPS6KA6 Protein (GST tag)

产品货号: TP08453

第三版

别名:PP90RSK4;RSK4

描述:Ribosomal protein S6 kinase alpha-6, also known as Ribosomal S6 kinase 4, 90 kDa ribosomal protein S6 kinase 6,RSK-4, RSK4 and RPS6KA6, is a protein which belongs to the protein kinase superfamily, AGC Ser/Thr protein kinase family and S6 kinase subfamily. RPS6KA6 contains one AGC-kinase C-terminal domain and two protein kinase domains. RPS6KA6 forms a complex with either ERK1 or ERK2 in quiescent cells. RPS6KA6 shows a high level of homology to three isolated members of the human RSK family. RSK2 is involved in Coffin-Lowry syndrome and nonspecific MRX. The localization of RPS6KA6 in the interval that is commonly deleted in mentally retarded males together with the high degree of amino acid identity with RSK2 suggests that RPS6KA6 plays a role in normal neuronal development. Further mutation analyses in males with X-linked mental retardation must prove that the gene of RPS6KA6 is indeed a novel MRX gene. RPS6KA6 is a serine/threonine kinase that may play a role in mediating the growth-factor and stress induced activation of the transcription factor CREB. RPS6KA6 is activated by multiple phosphorylations on threonine and serine residues.

配方:PBS

分子量:110 kDa

序列:Met 1-Leu 745

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.