

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human PCK2 Protein (His & GST Tag)

产品货号: TP08461

第三版

别名:PEPCK;PEPCK-M;PEPCK2

描述:PCK2 (Phosphoenolpyruvate Carboxykinase 2, Mitochondrial) is a Protein Coding gene. Diseases associated with PCK2 include Pepck 1 Deficiency and Phosphoenolpyruvate Carboxykinase Deficiency, Mitochondrial. Among its related pathways are Glycosaminoglycan metabolism and DNA Damage Response (only ATM dependent). Gene Ontology (GO) annotations related to this gene include GTP binding and purine nucleotide binding. An important paralog of this gene is PCK1.

配方:PBS

分子量:94.9 kDa

序列:Leu33-Met640

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.