

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human AK1 / Adenylate kinase 1 Protein (His tag)

产品货号: TP08478

第三版

别名:HTL-S-58j

描述:AK1 (Adenylate Kinase 1) is a Protein Coding gene. Diseases associated with AK1 include Adenylate Kinase Deficiency, Hemolytic Anemia Due To and Hemolytic Anemia. Among its related pathways are Metabolism and Pyrimidine metabolism (KEGG). Gene Ontology (GO) annotations related to this gene include kinase activity and adenylate kinase activity. An important paralog of this gene is AK5.

配方:PBS

分子量:23.5 kDa

序列:Met 1-Lys 194

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.