

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human CDKL2 Protein (His tag)

产品货号: TP08494

第三版

别名:KKIAMRE;P56

描述:Cyclin-dependent kinase-like 2 (Cdkl2) is a member of cdc2-related serine / threonine protein kinase family and it is found expressed in various brain regions, including the cerebral cortex, entorhinal cortex, hippocampus, amygdala, and dorsal thalamus. The high expression in these brain regions suggests that Cdkl2 may have functions in cognition and emotion processes. Cdkl2 accumulates primarily in the cytoplasm, with lower levels in the nucleus.

配方:PBS

分子量:58.3 kDa

序列:Met 1-His 493

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.