
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human MST4 Protein (GST tag)

产品货号: TP08509

第三版

别名:MASK;MST4

描述:MST4, also known as mammalian STE20-like protein kinase 4, is a novel member of the germinal center kinase subfamily of human Ste20-like kinases and is closely related to MST3. The 416 amino acid full-length MST4 contains a C-terminal regulatory domain and an N-terminal kinase domain, both of which are required for full activation of the kinase. MST4 is highly expressed in placenta, thymus, and peripheral blood leukocytes. MST4 specifically activates ERK but not JNK or p38 MAPK in transient transfected cells or in stable cell lines, and thus is biologically active in the activation of MEK/ERK pathway mediating cell growth and transformation. Further, MST4 kinase activity is stimulated significantly by epidermal growth factor receptor (EGFR) ligands, which are known to promote growth of certain cancer cells. Accordingly, MST4 have a potential role in signal transduction pathways involved in cancer progression. Three alternatively spliced isoform of MST4 have been isolated, and isoform 3 lacks an exon encoding kinase domain and may function as a dominant-negative regulator of the MST4 kinase.

配方:PBS**分子量:**73 kDa**序列:**Met 1-Pro 416**纯度:**> 95% by HPLC**浓度:****内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)**存储:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.