

14-3-3 ? Anti-Human Antibody (Clone J2E9)

Monoclonal Antibody

Cat.NO.: MA01016

3th Edition

Immunogen: Recombinant human 14-3-3 beta (1-246aa) purified from E. coli

Isotype: Mouse IgG2a heavy chain and kappa light chain

Clone: Anti-human 14-3-3 beta mAb, clone J2E9, is derived from hybridization of mouse FO myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human YWHAB protein.

Description:

Form: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

Concentration: 1 mg/ml

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage: The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application: ELISA, WB, ICC/IF