

**Alpha-Synuclein Antibody (Clone AT1E10)**

**Monoclonal Antibody**

**Cat.NO.: MA01043**

---

3th Edition

**Immunogen:** Recombinant human alpha-synuclein (1-140aa) purified from E. coli

**Isotype:** Mouse IgG1 heavy chain and kappa light chain

**Clone:** Anti-human alpha-synuclein mAb, clone AT1E10, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human alpha-synuclein protein.

**Description:**

**Form:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**Concentration:** 1 mg/ml

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:** The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:** ELISA, WB, ICC/IF