

**Amyloid Beta ( 11-28 ) ( 12B2 ) A? Anti-Human Mouse IgG MoAb**

**Monoclonal Antibody**

**Cat.NO.: MA01076**

---

3th Edition

**Immunogen:**Synthetic peptide for Human Amyloid [11-28]

**Isotype:**Mouse-Mouse hybridoma IgG1

**Clone:**12B2

**Description:**

**Form:**Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

**Concentration:**Please see datasheet.

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:**The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:**ELISA, WB, ICC/IF