

**Amyloid Beta ( 38-42 ) ( 44A3 ) A? Anti-Human Mouse IgG MoAb**

**Monoclonal Antibody**

**Cat.NO.: MA01080**

---

3th Edition

**Immunogen:** Synthetic peptide of human Amyloid Beta ( 38-42 ) (GVVIA).

**Isotype:** IgG2b

**Clone:** 44A3

**Description:**

**Form:** Lyophilized product in PBS containing 1 % BSA and 0.05 % NaN<sub>3</sub>.

**Concentration:** 1.0 mL deionized water will be added to the product, then its concentration comes to 50 µg/mL.

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:** The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:** ELISA, WB, ICC/IF