

**Amyloid Beta E22P ( 11A1 ) A? Anti-Human Mouse IgG MoAb**

**Monoclonal Antibody**

**Cat.NO.: MA01082**

---

3th Edition

**Immunogen:**Synthetic peptide of E22P- Amyloid Beta A? 10-35 part

**Isotype:**IgG1

**Clone:**11A1

**Description:**

**Form:**Lyophilized product from PBS containing 1 % BSA and 0.05 % NaN<sub>3</sub>

**Concentration:**Please see datasheet.

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:**The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:**ELISA, WB, ICC/IF