

**Angiotensinogen (104AT 601.2.80) Anti-Human Mouse IgG MoAb**

**Monoclonal Antibody**

**Cat.NO.: MA01084**

---

3th Edition

**Immunogen:** Recombinant protein of human angiotensinogen

**Isotype:** IgG1

**Clone:** 104AT 601.2.80

**Description:**

**Form:** Lyophilized product from PBS containing 1 % BSA and 0.05 % NaN<sub>3</sub>

**Concentration:** Please see datasheet.

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:** The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:** ELISA, WB, ICC/IF