

**Anti-Bcl-1/Cyclin D1 (5D4) Mouse IgG MoAb Biotin**

**Monoclonal Antibody**

**Cat.NO.: MA01109**

---

3th Edition

**Immunogen:**recombinant Human PRAD1/cyclin D1 (E. coli)

**Isotype:**IgG2a

**Clone:**5D4

**Description:**

**Form:**Lyophilized product from 1% BSA in PBS containing 0.05%NaN<sub>3</sub>

**Concentration:**100 ug/mL

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:**The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:**ELISA, WB, ICC/IF