

COX-2 (13H14) Anti-Human Mouse IgG MoAb

Monoclonal Antibody

Cat.NO.: MA01191

3th Edition

Immunogen: Synthetic peptide for a part of Human COX-2

Isotype: IgG1

Clone: 13H14

Description:

Form: Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

Concentration: 100 ug/mL

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage: The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application: ELISA, WB, ICC/IF