

**Anti-DPD (ADPYDMAB) Mouse IgG MoAb**

**Monoclonal Antibody**

**Cat.NO.: MA01221**

---

3th Edition

**Immunogen:**Human DPD protein

**Isotype:**Please see datasheet.

**Clone:**ADPYDMAB

**Description:**

**Form:**Lyophilized product in PBS containing 1 % BSA and 0.05 % NaN<sub>3</sub>

**Concentration:**Please see datasheet.

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:**The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:**ELISA, WB, ICC/IF