

EEF1A1 Antibody (Clone AT23C11)

Monoclonal Antibody

Cat.NO.: MA01228

3th Edition

Immunogen: Recombinant human EEF1A1 (1-462aa) purified from E.coli

Isotype: Mouse IgG2a, k

Clone: Anti-human EEFA1A mAb, clone AT23C11, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human EEF1A1 protein.

Description:

Form: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

Concentration: 1 mg/ml

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage: The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application: ELISA, WB, ICC/IF