

GAPDH Antibody (Clone AT4E8)

Monoclonal Antibody

Cat.NO.: MA01271

3th Edition

Immunogen:GAPDH from rabbit muscle

Isotype:Mouse IgG2b, k

Clone:Anti-human GAPDH mAb, clone AT4E8, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human GAPDH protein.

Description:

Form:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

Concentration:1 mg/ml

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage:The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application:ELISA, WB, ICC/IF