

GFAP (G-25-8-3) Anti-Bovine Mouse IgG MoAb

Monoclonal Antibody

Cat.NO.: MA01278

3th Edition

Immunogen:Bovine GFAP

Isotype:IgG2b

Clone:G-25-8-3

Description:

Form:Lyophilized product containing 0.05%NaN₃

Concentration:Please see datasheet.

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage:The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application:ELISA, WB, ICC/IF