

Glysyl-tRNA synthetase (GRS)

Monoclonal Antibody

Cat.NO.: MA01282

3th Edition

Immunogen: Full length native human recombinant GRS.

Isotype: IgG

Clone:

Description:

Form: Anti-Glysyl-tRNA synthetase antibody containing 0.05% sodium azide for preservative.

Concentration: Please see datasheet.

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage: The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application: ELISA, WB, ICC/IF