

**Glycyl-tRNA synthetase (GRS)**

**Monoclonal Antibody**

**Cat.NO.: MA01282**

---

3th Edition

**Immunogen:** Full length native human recombinant GRS.

**Isotype:** IgG

**Clone:**

**Description:**

**Form:** Anti-Glycyl-tRNA synthetase antibody containing 0.05% sodium azide for preservative.

**Concentration:** Please see datasheet.

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:** The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:** ELISA, WB, ICC/IF