

GRP94 Antibody (Clone AT94B9)

Monoclonal Antibody

Cat.NO.: MA01295

3th Edition

Immunogen: Recombinant human GRP94 (22-803aa) purified from E.coli.

Isotype: Mouse IgG2b heavy chain and kappa light chain.

Clone: Anti-human GRP94 mAb, clone AT94B9, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human GRP94 protein.

Description:

Form: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol

Concentration: 1 mg/ml

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage: The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application: ELISA, WB, ICC/IF