

Hepatic Sinusoidal Endothelial Cells (SE-1) Anti-Rat Mouse IgG MoAb

Monoclonal Antibody

Cat.NO.: MA01305

3th Edition

Immunogen:Rat Hepatic Sinusoidal Endothelial Cells

Isotype:IgG2a

Clone:SE-1

Description:

Form:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

Concentration:100 ug/mL

Storage: Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

Usage:The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

Application:ELISA, WB, ICC/IF