

**Anti-Human/Mouse/Rat ERK 1 Monoclonal Antibody**

**Monoclonal Antibody**

**Cat.NO.: MA02728**

---

3th Edition

**Immunogen:** Recombinant Protein of ERK1 of MAPK3

**Isotype:** IgG

**Clone:** 1C11

**Description:**

**Form:**

**Concentration:** 1mg/mL

**Storage:** Can be stored at +4°C. For long term storage, aliquot and store at -20°C. Avoid repeated freezing and thawing cycles.

**Usage:** The antibody has been tested by ELISA, Western blot analysis and ICC/IF to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results.

**Application:** ELISA, WB, ICC/IF