

c-Ret (Long Isoform) Anti-Human Rabbit IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01082

3th Edition

Description: The ret proto-oncogene products (proto-Ret protein) are expressed as 150kDa and 170kDa glycoproteins in neuroblastoma cells and as 150kDa and 190kDa glycoproteins in leukemia cells.

Antigen: Synthetic peptide for C-Terminal of Human c-Ret Long Isoform (ANWMLSPSAAKLMDTFDS)

Form: Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Long Isoform specific