

Anti-Endothelin-1 Rabbit IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01140

3th Edition

Description:Endothelins are 21 amino acid peptides that are produced ubiquitously by vascular endothelial cells, smooth muscle cells, and other cells in different organs. Endothelins are secreted as big-endothelins that are converted to active proteins by the endothelin-converting enzyme.

Antigen:Synthetic peptide of a part of ET-1 (A common parts of Bovine, Canine, Human, Mouse, Porcine and Rat)

Form:Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: