

EPRS (Glutamyl-prolyl tRNA synthetase)

Polyclonal Antibody

Cat.NO.: PA01146

3th Edition

Description: EPRS is a component of ARS multicomplex. EPRS is released and phosphorylated by IFN stimulation and forms a complex with GAIT.

Antigen: V573-K889 native human recombinant EPRS.

Form: Anti-glutamyl-prolyl-tRNA synthetase rabbit serum containing 0.05% sodium azide for preservative.

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: This antibody reacts specifically with 155kDa glutamyl prolyl tRNA synthetase protein.