

ERA (E. coli Ras-like protein) Anti-Human Rabbit IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01147

3th Edition

Description:ERA (Escherichia coli Ras-like protein) is an E. coli GTP binding protein that is essential for proliferation. A DNA database search suggests that homologous sequences with ERA exist in various organisms including human, mouse, Drosophila, Caenorhabditis elegans and Antirrhinum majus.

Antigen:Synthetic peptide of middle part of human ERA (E. coli Ras-like protein) ((C)AVKDPNTQSVGNPQ)

Form:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Confirmed by western blotting using transfectant cells