

**Fatty Acid Synthase Anti-Human Rabbit IgG Affinity Purify
Polyclonal Antibody**

Cat.NO.: PA01148

3th Edition

Description:Animal fatty acid synthase (FAS) is a homodimer protein, which synthesizes long-chain fatty acids and is rich in liver, brain, breast and lung. However, the precise cellular localization of FAS in human tissue has not been elucidated.

Antigen:Synthetic peptide for a part of Human Fatty Acid Synthase.

Form:Lyophilized product from 1 % BSA in PBS containing 0.05 % NaN₃

How to use:1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution:PBS (pH7.4) containing 1% BSA

Application:This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity:Human FAS specific. Cross reactivity with rat and mice are not tested