

**GRO/CINC-3 Anti-Rat Rabbit IgG**

**Polyclonal Antibody**

**Cat.NO.: PA01187**

---

3th Edition

**Description:** Cytokine-induced neutrophil chemo attractant 1 (CINC-1) was originally purified from media conditioned by IL-1  $\alpha$  stimulated rat kidney epithelioid cells (NRK-52E).

**Antigen:** Recombinant Rat GRO/CINC-3.

**Form:** Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1  $\mu$ g/ml.

**Specificity:** Rat GRO/CINC-3 specific. Non-cross react with Rat GRO/CINC-1, Rat GRO/CINC-2 $\alpha$  and Rat GRO/CINC-2 $\beta$ .