

HGF? (H56) Anti-Rat Rabbit IgG Affinity Purify

Polyclonal Antibody

Cat.NO.: PA01247

3th Edition

Description: Hepatocyte Growth factor (HGF) was discovered as a mitogen for hepatocytes. HGF was subsequently found to be identical to the scatter factor, which destroys epithelial cell adhesion and promotes cell movement.

Antigen: Synthetic peptide for N-Terminal of rat HGF

Form: Lyophilized product from 1% BSA in PBS containing 0.05%NaN₃

How to use: 1.0 ml distilled water will be added to the product

Stability: Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

Dilution: PBS (pH7.4) containing 1% BSA

Application: This antibody can be used for western blotting in concentration of 1?5?g/ml.

Specificity: Rat HGF-? specific. Non-cross react with Human HGF and Rat HGF?.