

**HGF? (H56) Anti-Rat Rabbit IgG Affinity Purify**

**Polyclonal Antibody**

**Cat.NO.: PA01247**

---

3th Edition

**Description:** Hepatocyte Growth factor (HGF) was discovered as a mitogen for hepatocytes. HGF was subsequently found to be identical to the scatter factor, which destroys epithelial cell adhesion and promotes cell movement.

**Antigen:** Synthetic peptide for N-Terminal of rat HGF

**Form:** Lyophilized product from 1% BSA in PBS containing 0.05%NaN<sub>3</sub>

**How to use:** 1.0 ml distilled water will be added to the product

**Stability:** Lyophilized product, 5 years at 2 – 8°C; Solution, 2 years at –20°C

**Dilution:** PBS (pH7.4) containing 1% BSA

**Application:** This antibody can be used for western blotting in concentration of 1?5?g/ml.

**Specificity:** Rat HGF-? specific. Non-cross react with Human HGF and Rat HGF?.